

GVL

Banc gamma vertical

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GVL

Vertical gammadensitometer bench

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Apparatus for the investigation of materials of various types by measurement of their density using the absorption of gamma rays. It can be used to check the homogeneity of a specimen, track changes in density during compaction, etc.; it can also be used to measure the densities of materials taken in containers (mud, clays, silts, etc.).

The apparatus includes:

- a source and a detector on the same vertical axis, rigidly connected together and capable of moving 345 mm horizontally,
- a carriage, holding the source and the detector, that can move 1,000 mm longitudinally;
- a frame holding the moving parts and the specimens holder (mobile plate for prismatic specimens or motor-driven rollers with longitudinal axes for cylindrical specimens);
- systems for measuring the thickness of material examined, metering the radiation, and controlling the movements.

It allows automatic examination of up to 999 points on prismatic specimens (point-by-point mode), or continuously along a chord or helix, on cylindrical specimens (continuous mode). A PC-type microcomputer controls the automatic functions and stores the parameters; a pre-processing program can be used to generate an ASCII file of the apparent densities and void contents.

Sample size:

- prism: $h < 400$ mm; mass < 150 kg
- cylinder: $50 < \varnothing < 200$ mm ; $l < 1,000$ mm; mass < 75 kg

Equipment size: $w = 1.95$ m; $d = 1.47$ m to 1.60 m; $h = 1.97$ m; Mass: 600 kg

Counting system size: $w = 0.25$ m; $d = 0.40$ m; $h = 0.15$ m; Mass: 7 kg